

Validate the interface between our L.E.S.S. lighting systems and your microscope by using the below table.

Olympus											
Microscope Series	Microscope Type	Objective Olympus #	Objective Type	Mag	Diameter in mm	Working Distance in mm	Brightfield L.E.S.S. VisiBright	Darkfield L.E.S.S. VisiDark	Ring Adapter L.E.S.S. AC-RA-BF-58-80 or AC-RA-BF-50-58-C	Combo Post L.E.S.S. AC-P-M	
SZ2	SZ 51 SZ 61	-	-	1x	54	110	VisiBright	VisiDark	AC-RA-BF-50-58-C	AC-P-M	
		110AL0,25x	Achro	0,25x	54	400	-	VisiDark	-	-	
		110AL0,5x	Achro	0,5x	54	200	-	VisiDark	-	-	
		110AL0,62x	Achro	0,62x	54	160	-	VisiDark	-	-	
		110AL0,75x	Achro	0,75x	54	130	VisiBright	VisiDark	AC-RA-BF-50-58-C	AC-P-M	
		110AL1,5x	Achro	1,5x	54	61	VisiBright	VisiDark	AC-RA-BF-50-58-C	AC-P-M	
		110AL2x	Achro	2x	54	38	-	VisiDark	AC-RA-BF-50-58-C	AC-P-M	
SZX2	SZX7 SZX10	DFPL0,5x-4	Plan	0,5x	60	171	-	VisiDark	-	-	
		DFPL0,75x-4	Plan	0,75x	60	116	VisiBright	VisiDark	AC-RA-BF-58-80	AC-P-M	
		SZX-ACH1x	Achro	1x	60	90	VisiBright	VisiDark	AC-RA-BF-58-80	AC-P-M	
		SZX-ACH1,25x-2	Achro	1,25x	60	68	VisiBright	VisiDark	AC-RA-BF-58-80	AC-P-M	
		DFPLAPO1,25x-4	Planapo	1,25x	60	60	VisiBright	VisiDark	AC-RA-BF-58-80	AC-P-M	
		DFPLAPO1x-4	Planapo	1x	60	81	VisiBright	VisiDark	AC-RA-BF-58-80	AC-P-M	
		DFPL1,5x4	Plan	1,5x	60	45,5	-	VisiDark	AC-RA-BF-58-80	AC-P-M	
		DFPL2x4	Plan	2x	60	33,5	-	VisiDark	AC-RA-BF-58-80	AC-P-M	
		SZX16	SDFPLFL0,3X	Plan	0,3x	60	141	VisiBright	VisiDark	AC-RA-BF-58-80	AC-P-M
	SDFPLAPO0,5XPF		Planapo	0,5x	60	70,5	VisiBright	VisiDark	AC-RA-BF-58-80	AC-P-M	
	SDFPLAPO0,8X		Planapo	0,8x	60	81	VisiBright	VisiDark	AC-RA-BF-58-80	AC-P-M	
	SDFPLAPO1XPF		Planapo	1x	60	60	VisiBright	VisiDark	AC-RA-BF-58-80	AC-P-M	
	SDFPLAPO1,6XPF		Planapo	1,6x	68,5	30	-	VisiDark	AC-RA-BF-58-80	AC-P-M	
	SDFPLAPO2XPFC		Planapo	2x	74	20	-	VisiDark	AC-RA-BF-58-80	AC-P-M	
SZX	SZX12	DFPLFL0,3x	Plan	0,3x	60	130	VisiBright	VisiDark	AC-RA-BF-58-80	AC-P-M	
		DFPLFL0,45x	Plan	0,45x	60	198	-	VisiDark	-	-	
		DFPLFL0,5x-PF	Plan	0,5x	60	70	VisiBright	VisiDark	AC-RA-BF-58-80	AC-P-M	
		DFPLAPO1x-PF	Planapo	1x	60	74	VisiBright	VisiDark	AC-RA-BF-58-80	AC-P-M	
		DFPLAPO1,2x-PF2	Planapo	1,2	60	60	VisiBright	VisiDark	AC-RA-BF-58-80	AC-P-M	
		DFPLFL1,6x-PF	Plan	1,6	60	34	-	VisiDark	AC-RA-BF-58-80	AC-P-M	
		SZX-AL20	Extension	2,5x	-	10	-	VisiDark	-	-	
			SZX9	DFPL0,5x	Plan	0,5x	60	198	-	VisiDark	-
	DFPL0,75x	Plan		0,75x	60	113	VisiBright	VisiDark	AC-RA-BF-58-80	AC-P-M	
	DFPLAPO1x-2	Planapo		1x	60	87,5	VisiBright	VisiDark	AC-RA-BF-58-80	AC-P-M	
	SZX-ACH1x	Achro		1x	60	90	VisiBright	VisiDark	AC-RA-BF-58-80	AC-P-M	
	DFPL1,5x	Plan		1,5x	60	53	VisiBright	VisiDark	AC-RA-BF-58-80	AC-P-M	
	DFPL2x-3	Plan		2x	60	34	-	VisiDark	AC-RA-BF-58-80	AC-P-M	